



# Solutions for Family Law

## Prove Sobriety—Maintain Custody and Visitation Rights



In divorce, custody, or visitation cases, when a parent's alcohol issues are called into question, the ability to care for the children in the home becomes a concern. Judges, family law attorneys, alternate dispute arbitrators, or social service agencies can regularly require parents to take sobriety tests. A parent involved in a custody or visitation arrangement may also volunteer to validate sobriety with BI SL3®. Using the right alcohol detection tool, results in sound custody and visitation decisions. The SL3 handheld breathalyzer makes it easy to monitor sobriety, while driving compliance and accountability.



### Benefits of Mobile Alcohol Monitoring

The SL3 is a user-friendly, lightweight alcohol monitoring device that discreetly prompts the user to test. This is especially convenient for individuals who may be at work or in a social setting. The SL3 uses the highest quality fuel-cell sensor to achieve an accuracy level within .005 BrAC. Plus, real-time delivery of sobriety results offers immediate peace of mind for families in custody or visitation arrangements.

Receiving urine analysis results can take days, and negative results do not necessarily mean abstinence, as alcohol can quickly be eliminated from the body. With SL3, multiple tests per day can be conducted at scheduled or random times per the custody or visitation agreement. Notifications are sent to the client's mobile device, prompting him or her to submit a test. Test results are then conveniently sent to the client's mobile phone or email in real-time.



## Key Features

- Discreet handheld device
- Test reminders via text message
- .02 BrAC or greater prompts a re-test
- Customizable test scheduling
- Real-time text or email alerts for violations
- Automated daily, weekly, or monthly reports
- Biometric facial comparison technology confirms client identity
- Professional grade fuel cell
- Backlit LCD screen for easy reading of prompts and test taking
- Stores 70 tests

## How it Works

- Client enrolls, and testing schedule and alerts are set up
- Automated text message reminders are sent to the client's mobile device based on schedule
- Client tests within 45 minutes of being prompted
- Real-time results and alerts on missed or positive results are sent to all designated recipients



For more information, call 800.876.1105 or visit us at [bi.com](http://bi.com)